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Ron Vale's laboratory at UCSF and is funded by an NIH grant R01-EB007187 from the National Institute of Biomedical Imaging and Bioengineering (NBIB).

Arthur Edelstein, Nenad Amodaj, Karl Hoover, Ron Vale, and Nico Stuurman

Department of Cellular and Molecular Pharmacology, University of California

San Francisco, San Francisco, California


100x Imaging, San Francisco, California

Howard Hughes Medical Institute, Chevy Chase, Maryland

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




L'interface graphique

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
Installed Devices:


Name	Adapter/Library	Description	Status
Core	MMCore/Default	Core controller	Default
DHub	DHub/DemoCamera	DHub	OK
Camera	DCam/DemoCamera	Demo camera	OK
Dichroic	DWheel/DemoCam...	Demo filter wheel	OK
Emission	DWheel/DemoCam...	Demo filter wheel	OK
Excitation	DWheel/DemoCam...	Demo filter wheel	OK
Objective	DObjective/DemoC...	Demo objective turret	OK
Z	DStage/DemoCamera	Demo stage	OK
Path	DLightPath/DemoC...	Demo light path	OK
XY	DXYStage/DemoCa...	Demo XY stage	OK
Shutter	DShutter/DemoCa...	Demo shutter	OK
Autofocus	DAutoFocus/Demo...	Demo auto focus	OK

Available Devices: list by vendor compact view

- AAAOTF
- ABSCamera
- AOTF
- ASIFW1000
- ASISStage
- ASIwptr
- Aladdin
- Andor
- AndorLaserCombiner
- Aquinas
- Arduino
- BaumerOptron
- CARVII
- CSUX
- Cobolt
- CelbescentCuba

La fenêtre de script en beanshell

Exemple: Media:acqLC.bsh - complicated MD acquisition including visiting multiple positions as defined in the position list, acquiring time-lapse data in individual channels.

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MMStudio Version 1.4.6 20110907

Configuration file:

C:\Program Files\Micro-Manager1.4-20110731\MMConfig_demo...

This software is distributed free of charge in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of merchantability or fitness for a particular purpose. In no event shall the copyright owner or contributors be liable for any direct, indirect, incidental, special, exemplary, or consequential damages.




Copyright University of California San Francisco, 2007, 2008, 2009, 2010. All rights reserved.

Micro-Manager was initially funded by grants from the Sandler Foundation and is now supported by a grant from the NIH.

If you have found this software useful, please cite Micro-Manager in your publications.

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SutterLambda

Summary: Interfaces to Sutter Lambda 10 controllers, Lambda SC controller, and DG4

Authors: Nenad Amodaj

Maintainers: Arthur Edelstein

License: BSD

Platforms: Windows, Mac OS X and Linux (only serial port, no parallel port or USB)

Devices: Lambda SC

Lambda 10-B

Lambda 10-2

Lambda 10-3

Lambda DG-4

Contents

1 Device Configuration

1.1 Serial port issues

1.2 Wheel speed

2 Troubleshooting

2.1 Wheel re-homing during imaging


Device Configuration

Serial port issues

Sutter devices use a straight through (i.e., not null modem) DB9 serial cable. The default baud rates are 9600 for serial ports and 128800 for USB ports (which use a build-in USB-serial converter)

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State devices

Roue_Filtres

State	Label
0	FITC
1	TexasRed
2	Cy5
3	Vide
4	Filter-4
5	Filter-5
6	Filter-6
7	Filter-7
8	Filter-8
9	Filter-9

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Configuration settings Save

Group	Preset
Camera	Default
Channel	Rhodamine
Gain	5 < >
LightPath	Eyepiece
Objective	40X

Group: + - Edit Preset: + - Edit

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Configuration

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Group Editor

Here you can specify the properties included in a configuration group.

☐ Show read-only properties

Group name:

Show

- ☒ cameras
- ☒ shutters
- ☒ stages
- ☒ wheels, turrets, etc.
- ☒ other devices

Property Name | **Use in Group?** | **Current Property Value**

Camera-TestProperty1	<input type="checkbox"/>	0
Camera-TestProperty5	<input type="checkbox"/>	0
Camera-TestProperty6	<input type="checkbox"/>	0
Camera-TransposeCorrection	<input type="checkbox"/>	0
Camera-TransposeMirrorX	<input type="checkbox"/>	0
Camera-TransposeMirrorY	<input type="checkbox"/>	0
Camera-TransposeXY	<input type="checkbox"/>	0
Camera-TriggerDevice	<input type="checkbox"/>	0
Camera-UseExposureSequences	<input type="checkbox"/>	No
Roue_Filtres-Busy	<input type="checkbox"/>	0
Roue_Filtres-ClosedPosition	<input type="checkbox"/>	0
Roue_Filtres-Delay_ms	<input type="checkbox"/>	0
Roue_Filtres-Label	<input checked="" type="checkbox"/>	TexasRed
Roue_Filtres-Speed	<input type="checkbox"/>	3
Roue_Filtres-State	<input type="checkbox"/>	1
Core-AutoFocus	<input type="checkbox"/>	1
Core-AutoShutter	<input type="checkbox"/>	1
Core-Camera	<input type="checkbox"/>	Camera
Core-ChannelGroup	<input type="checkbox"/>	
Core-Focus	<input type="checkbox"/>	
Core-Galvo	<input type="checkbox"/>	
Core-ImageProcessor	<input type="checkbox"/>	
Core-Initialize	<input type="checkbox"/>	1
Core-SLM	<input type="checkbox"/>	
Core-Shutter	<input type="checkbox"/>	
Core-TimeoutMs	<input type="checkbox"/>	5000

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Device -- loads device

Label -- attaches a label to the specified device position

Equipment -- defines various equipment attributes used for image annotation

Property -- sets a specific device property

ConfigGroup -- defines a single entry in the configuration group

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- Exemple

Generated by Configurator on Mon Jun 01 11:56:47
CEST 2015

Reset
Property,Core,Initialize,0

Devices
Device,COM1,SerialManager,COM1
Device,Roue_Filtres,SutterLambda,Wheel-A
Device,DHub,DemoCamera,DHub
Device,DCam,DemoCamera,Dcam

Pre-init settings for devices
Property,Roue_Filtres,Port,COM1
Property,DCam,MaximumExposureMs,10000.0000

Pre-init settings for COM ports
Property,COM1,AnswerTimeout,500.0000
Property,COM1,BaudRate,9600
Property,COM1,DelayBetweenCharsMs,0.0000
Property,COM1,Handshaking,Off
Property,COM1,Parity,None
Property,COM1,StopBits,1
Property,COM1,Verbose,1

Hub (parent) references
Parent,DCam,Dhub
Initialize
Property,Core,Initialize,1

Delays
Roles
Property,Core,Camera,DCam
Property,Core,AutoShutter,1

Camera-synchronized devices
Labels
Roue_Filtres
Label,Roue_Filtres,3,Vide
Label,Roue_Filtres,2,Cy5
Label,Roue_Filtres,1,TexasRed
Label,Roue_Filtres,0,FITC

Configuration presets
Group: Roue
Preset: TexasRed
ConfigGroup,Roue,TexasRed,Roue_Filtres,Label,TexasRed

Les outils plus avancés

- Le MDA

Multi-dimensional Acquisition

☒ **Time points**

Number: 10

Interval: 0 ms

☒ **Multiple positions (XY)**

Edit position list...

☒ **Z-stacks (slices)**

Z-start [um]: 0 Set

Z-end [um]: 0 Set

Z-step [um]: 1

relative Z

☐ Keep shutter open

☒ **Channels**

Channel group: Channel ☐ Keep shutter open

Use?	Configu...	Exposure	Z-offset	Z-stack	Skip Fr.	Color
<input checked="" type="checkbox"/>	Cy5	10	0	<input checked="" type="checkbox"/>	0	
<input type="checkbox"/>	DAPI	10	0	<input checked="" type="checkbox"/>	0	

New Remove Up Down

☒ **Save images**

Directory root: C:\Users\ ...

Name prefix: Untitled

Acquisition Comments

Summary

Number of time points: 10
Number of positions: 1
Number of slices: 1
Number of channels: 1
Total images: 10
Total memory: 5 MB
Duration: 0h 0m 0s
Order: Time, Position, Channel, Slice

Advanced

Close

Acquire!

Stop

Load...

Save as...